

Lambda CE6 Induction Kit

INSTRUCTION MANUAL

Catalog #235200

Revision A

For In Vitro Use Only

235200-12

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Agilent. Agilent shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

ORDERING INFORMATION AND TECHNICAL SERVICES

United States and Canada

Agilent Technologies

Stratagene Products Division

11011 North Torrey Pines Road

La Jolla, CA 92037

Telephone (858) 373-6300

Order Toll Free (800) 424-5444

Technical Services (800) 894-1304

Internet techservices@agilent.com

World Wide Web www.stratagene.com

Europe

Location	Telephone	Fax	Technical Services
Austria	0800 292 499	0800 292 496	0800 292 498
Belgium	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 15775	0800 15740	0800 15720
France	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 919 288	0800 919 287	0800 919 289
Germany	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 182 8232	0800 182 8231	0800 182 8234
Netherlands	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 023 0446	+31 (0)20 312 5700	0800 023 0448
Switzerland	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 563 080	0800 563 082	0800 563 081
United Kingdom	00800 7000 7000	00800 7001 7001	00800 7400 7400
	0800 917 3282	0800 917 3283	0800 917 3281

All Other Countries

Please contact your local distributor. A complete list of distributors is available at www.stratagene.com.

Lambda CE6 Induction Kit

CONTENTS

Materials Provided.....	1
Storage Conditions.....	1
Additional Materials Required	1
Academic and Nonprofit Laboratory Assurance Letter.....	2
Introduction.....	3
Host and Expression Strains	3
LE392 Host Cells	3
BL21 Competent Cells	3
Amplifying the Lambda CE6 Bacteriophage.....	4
Plate Lysate Amplification Method.....	4
Titering the Lambda CE6 Bacteriophage.....	5
Transformation Guidelines	6
Storage Conditions	6
Aliquoting Cells	6
Use of 14-ml BD Falcon Polypropylene Tubes.....	6
Use of β -Mercaptoethanol.....	6
Quantity of DNA Added	6
Length of the Heat Pulse	6
Transformation Protocol.....	7
Induction of Target Protein by Infection with Bacteriophage Lambda CE6	9
Troubleshooting	10
Preparation of Media and Reagents.....	11
References	12
MSDS Information.....	12

Lambda CE6 Induction Kit

MATERIALS PROVIDED

Materials provided	Quantity
Bacteriophage lambda CE6 (titer is $\geq 5.0 \times 10^9$ pfu/ml)	1.0 ml from high-titer lambda lysate stock
LE392 host cells (bacterial glycerol stock)	1.0 ml
BL21 competent cells	Five 0.2-ml tubes
pUC18 control plasmid ^a	10 μ l
β -Mercaptoethanol (1.42 M)	25 μ l

^a The pUC18 control plasmid is shipped at a concentration of 0.1 ng/ μ l in TE buffer (see *Preparation of Media and Reagents*).

STORAGE CONDITIONS

All Materials: -80°C

ADDITIONAL MATERIALS REQUIRED

14-ml BD Falcon polypropylene round-bottom tubes (BD Biosciences Catalog #352059)

Revision A

© Agilent Technologies, Inc. 2008.

ACADEMIC AND NONPROFIT LABORATORY ASSURANCE LETTER

The T7 expression system is based on technology developed at Brookhaven National Laboratory under contract with the U.S. Department of Energy and is protected by U.S. patents assigned to Brookhaven Science Associates (BSA). BSA will grant a nonexclusive license for use of this technology, including the enclosed materials, based on the following assurances:

1. These materials are to be used for noncommercial research purposes only. A separate license is required for any commercial use, including the use of these materials for research purposes or production purposes by any commercial entity. Information about commercial licenses may be obtained from the Office of Intellectual Property and Industrial Partnerships, Brookhaven National Laboratory, Bldg. 475D, Upton, New York, 11973 [telephone (631) 344-7134].
2. No materials that contain the cloned copy of T7 gene 1, the gene for T7 RNA polymerase, may be distributed further to third parties outside of your laboratory, unless the recipient receives a copy of this license and agrees to be bound by its terms. This limitation applies to strains BL21, BL21(DE3), and BL21(DE3)pLysS included in this kit and any derivatives you may make of them.

You may refuse this license by returning the enclosed materials unused. By keeping or using the enclosed materials, you agree to be bound by the terms of this license.

Commercial Entities Outside of the US

The T7 expression system is based on technology developed at Brookhaven National Laboratory under contract with the U.S. Department of Energy and is protected by U.S. Patents assigned to Brookhaven Science Associates (BSA). To protect its patent properties BSA requires commercial entities doing business in the United States, its Territories or Possessions to obtain a license to practice the technology. This applies for in-house research use of the T7 system as well as commercial manufacturing using the system. Commercial entities outside the U.S. that are doing business in the U.S., must also obtain a license in advance of purchasing T7 products. Commercial entities outside the U.S. that are using the T7 system solely for in-house research need not obtain a license if they do no business in the United States. However all customers, whether in the U.S. or outside the U.S. must agree to the terms and conditions in the Assurance Letter which accompanies the T7 products. Specifically, no materials that contain the cloned copy of T7 gene 1, the gene for T7 RNA polymerase, may be distributed further to third parties outside of your laboratory, unless the recipient receives a copy of the assurance letter and agrees to be bound by its terms. This limitation applies to strains BL21, BL21(DE3), and BL21(DE3)pLysS included in this kit and any derivatives you may make of them.

To obtain information about licensing, please contact the Office of Intellectual Property and Industrial Partnerships, Brookhaven National Laboratory, Building 475D, Upton, NY 11973 [telephone: 631-344-7134; Fax: 631-344-3729].

INTRODUCTION

The Lambda CE6 Induction Kit* is used to deliver T7 RNA polymerase into BL21 competent cells to induce expression of toxic proteins from T7 promoter-driven vectors, such as the pCAL and pET vectors. The Lambda CE6 Induction Kit allows tight control of protein expression, which is important when expressing toxic proteins since no T7 RNA polymerase is present in the cells until the desired time of induction. Bacteriophage lambda CE6 expresses T7 RNA polymerase from the lambda p_L and p_I promoters and carries the *Sam7* lysis mutations. T7 polymerase drives the transcription of the gene on the expression plasmid downstream of the T7 promoter. Bacteriophage lambda CE6 allows effective expression of target genes in BL21 cells and presumably other nonrestricting hosts which adsorb lambda. Lambda bacteriophage CE6 is propagated in the LE392 host strain,¹ which suppresses the *Sam7* mutation and consequently undergoes lysis.

HOST AND EXPRESSION STRAINS

LE392 Host Cells

e14⁻(mcrA) hsdR514 supE44 supF58 lacY1 or Δ(lacIZY)6 galK2 galT22 metB1 trpR55

Streak bacterial glycerol stock on either LB or NZY agar plates (see *Preparation of Media and Reagents*).

BL21 Competent Cells

E. coli B F⁻ dcm ompT hsdS(r_B⁻ m_B⁻) gal

See reference 2.

* U.S. Patent No. 4,952,496.

AMPLIFYING THE LAMBDA CE6 BACTERIOPHAGE

It is important to titer the lambda CE6 bacteriophage prior to each use. Expect titers of approximately 5×10^9 pfu/ml. If the titer drops over time, prepare a fresh high-titer stock of the lambda bacteriophage as outlined in the amplification method below. For further details about lambda phage titering or amplification methods, see reference 3.

Plate Lysate Amplification Method

1. Transfer a single colony of LE392 cells into 5 ml of NZY broth with maltose.[§] Incubate with shaking at 37°C until growth reaches an OD₆₀₀ of 1.0.
2. Centrifuge the overnight culture at 1700–2000 × g for 15 minutes at 4°C. Discard the supernatant. Resuspend the cell pellet in 10 mM MgSO₄ to a final OD₆₀₀ of 0.5.
3. Combine 250 µl of cells (at OD₆₀₀ = 0.5) with 1×10^6 pfu of CE6 in 14-ml BD Falcon polypropylene round-bottom tubes in triplicate. Incubate the tubes at 37°C for 15 minutes without shaking to allow the phage to attach to the cells.
4. Add 3 ml of NZY top agar,[§] melted and cooled to ~48°C, to each cell suspension and plate on separate dry, prewarmed 100-mm agarose plates.[§] Allow the plates to set for 10 minutes. Invert the plates and incubate overnight at 37°C.
5. Overlay each plate with 5 ml of SM buffer.[§] Rock the plates for 4 hours at room temperature to allow the phage to diffuse into the SM buffer.
6. Remove the SM buffer (which contains the lambda CE6) from each plate and pool the volumes in a 50-ml conical tube.
7. Add 500 µl of chloroform. Mix well and incubate for 15 minutes at room temperature.
8. Centrifuge the SM buffer at 1700–2000 × g for 15 minutes at 4°C to pellet the cell debris and chloroform.
9. Remove the phage-containing supernatant and determine the titer of the solution (see *Titering the Lambda CE6 Bacteriophage*). Store the lambda CE6 stock at 4°C.

[§] See *Preparation of Media and Reagents*.

TITERING THE LAMBDA CE6 BACTERIOPHAGE

1. Transfer a single colony of LE392 cells into 50 ml of NZY broth with maltose.
2. Incubate with shaking at 37°C for 4–6 hours (do not grow past an OD₆₀₀ of 1.0). Alternatively, grow overnight at 30°C, shaking at 200 rpm.

Note *The lower temperature keeps the bacteria from overgrowing, thus reducing the number of nonviable cells. Phage can adhere to nonviable cells resulting in a decreased titer.*

3. Pellet the bacteria at 1000 × g for 10 minutes.
4. Gently resuspend the cell pellet in sterile 10 mM MgSO₄ and dilute to an OD₆₀₀ of 0.5.

Note *The bacteria should be used immediately following dilution.*

5. Dilute aliquots of the lambda CE6 bacteriophage stock in SM buffer by the following amounts: 1:100, 1:1000, 1:10,000. Add 1 µl of each dilution to 200 µl of LE392 host cells at an OD₆₀₀ of 0.5.
6. Incubate the phage and the bacteria at 37°C for 15 minutes to allow the phage to attach to the cells.
7. Add 3 ml of NZY top agar, melted and cooled to ~48°C, and plate immediately onto dry, prewarmed agarose plates. Allow the plates to set for 10 minutes. Invert the plates and incubate at 37°C.
8. Plaques should be visible after 6–8 hours. Count the plaques and determine the titer in plaque-forming units per milliliter (pfu/ml).

TRANSFORMATION GUIDELINES

Important *To achieve optimal transformation efficiency, please read the guidelines outlined in the following sections before proceeding with the Transformation Protocol.*

Storage Conditions

The competent cells are very sensitive to slight variations in temperature. Storing the competent cells at the bottom of a -80°C freezer directly from the dry ice shipping container is required in order to prevent a loss of transformation efficiency. Transferring tubes from one freezer to another may result in a loss of efficiency. The transformation efficiency of the competent cells is guaranteed when the competent cells are used according to the specifications outlined in this instruction manual.

Aliquoting Cells

Store the competent cells on ice at all times while aliquoting. It is essential to place the BD Falcon polypropylene tubes on ice before the competent cells are thawed and to aliquot the competent cells directly into the prechilled BD Falcon polypropylene tubes. It is also important to use 100 μl of competent cells/transformation. Using an inadequate volume of competent cells results in lower transformation efficiencies.

Use of 14-ml BD Falcon Polypropylene Tubes

It is important to use 14-ml BD Falcon polypropylene round-bottom tubes for the *Transformation Protocol*, because other tubes may be degraded by β -mercaptoethanol. Additionally, the critical incubation period during heat-pulsing is calculated for the thickness and shape (i.e., round bottom) of the BD Falcon polypropylene tubes.

Use of β -Mercaptoethanol

β -Mercaptoethanol increases transformation efficiencies two- to threefold. Use 1.7 μl of β -mercaptoethanol provided with this kit or a fresh 1:10 dilution of a 14.2 M stock solution per 100 μl of cells.

Quantity of DNA Added

Greatest efficiencies (i.e., transformants per microgram of DNA) are observed when adding 1 μl of ligated DNA at a concentration of 0.1 ng/ μl per 100 μl of competent cells. Although the overall transformation efficiency may be lower, a greater number of colonies will be obtained when transforming up to 50 ng.

Length of the Heat Pulse

Optimal transformation efficiencies are observed when transformation reactions are heat-pulsed for 45–50 seconds. Transformation efficiencies decrease sharply when heat-pulsed for <45 seconds or for >60 seconds.

TRANSFORMATION PROTOCOL

1. Thaw the BL21 competent cells on ice.

Note *Store the competent cells **on ice at all times** while aliquoting. It is essential that the BD Falcon polypropylene tubes are placed on ice before the competent cells are thawed and that 100 μ l of competent cells are aliquoted directly into each **prechilled** polypropylene tube. Pipet the remaining competent cells into 100- μ l aliquots and freeze the aliquots at -80°C . Do not pass the frozen competent cells through more than one freeze-thaw cycle.*

2. Gently mix the competent cells. Aliquot 100 μ l of the competent cells into the appropriate number of prechilled 14-ml BD Falcon polypropylene tubes.
3. Add 1.7 μ l of β -mercaptoethanol provided with this kit or a fresh 1:10 dilution of a 14.2 M stock solution of β -mercaptoethanol, diluted in distilled water (dH_2O), to each polypropylene tube containing the competent cells to a final concentration of 25 mM and swirl gently.
4. Incubate the reactions on ice for 10 minutes, swirling gently every 2 minutes.
5. Add 1–50 ng of ligated DNA to each transformation reaction and swirl gently. For the control transformation reaction, add 1 μ l of the pUC18 control plasmid to a separate 100- μ l aliquot of the competent cells and swirl gently.
6. Incubate the reactions on ice for 30 minutes.
7. Heat-pulse each transformation reaction in a 42°C water bath for 45 seconds. **The duration of the heat pulse is critical for optimal transformation efficiencies.**
8. Incubate the reactions on ice for 2 minutes.
9. Add 0.9 ml of SOC medium (see *Preparation of Media and Reagents*) to each transformation reaction and incubate the reactions at 37°C for 1 hour with shaking at 225–250 rpm.

10. Using a sterile spreader, plate ≤ 200 μl of each transformation reaction directly onto separate LB agar plates containing the appropriate antibiotic. If plating ≥ 100 μl , the cells can be spread directly onto the plates. If plating < 100 μl of the transformation reaction, plate into a 200- μl pool of SOC medium. When spreading bacteria onto the plate, tilt and tap the spreader to remove the last drop of cells.

Use a sterile spreader to plate 200 μl of the control transformation reaction containing the pUC18 plasmid directly onto an LB–ampicillin agar plate (see *Preparation of Media and Reagents*). The expected results for the control transformation reaction are as follows:

Host strain	Transformation plated	Expected colonies	Efficiency (cfu/ μg of pUC18 DNA)
BL21 strain	200 μl	> 20	$\geq 1 \times 10^6$

INDUCTION OF TARGET PROTEIN BY INFECTION WITH BACTERIOPHAGE LAMBDA CE6

Note *This protocol is designed for inductions in 50-ml culture volumes. If inductions of larger volumes of culture are desired, it will be necessary to increase the volume of the overnight culture in step 1. The increased volume of overnight culture is necessary to achieve the required cell density ($A_{600} \leq 1$) in the larger volume of broth the following day.*

1. Inoculate 5 ml of NZY broth with maltose containing the antibiotic required to maintain the expression plasmid with a single colony of BL21 cells harboring the expression plasmid. Shake overnight at 37°C at 200–250 rpm.

Note *If the competent cells contain a pACYC-based plasmid (e.g., any BL21-CodonPlus strain), the overnight culture must include chloramphenicol at a final concentration of 50 µg/ml in addition to the antibiotic required to maintain the expression plasmid.*

2. In the morning, centrifuge 1.0 ml of the overnight culture, resuspend the cells in 1.0 ml of fresh NZY broth with maltose, and pipet the resuspended cells into a flask containing 50 ml of fresh NZY broth with maltose.
3. Record the A_{600} of the diluted culture. It should be ≤ 0.1 . If the A_{600} is >0.1 , dilute the culture with fresh NZY broth with maltose to an A_{600} of ≤ 0.1 . If the A_{600} is <0.1 , the time required to reach an A_{600} of 0.3 (in step 4) will be extended.
4. Grow the culture to an A_{600} of 0.3 and add glucose to a final concentration of 4 mg/ml (e.g., 1.0 ml of a 20% glucose solution to the 50-ml culture).
5. Grow the culture to an A_{600} of 0.6–1.0 and add MgSO_4 to a final concentration of 10 mM (e.g., 500 µl of 1.0 M MgSO_4 to the 50-ml culture).

6. Remove a portion of the culture to serve as the uninduced control and infect the rest with bacteriophage lambda CE6 at a multiplicity of infection (MOI) of 5–10 particles per cell. (To optimize induction, cultures may be split into 3 or 4 aliquots and infected with varying dilutions of bacteriophage lambda CE6. The subsequent induction can be monitored by SDS-PAGE or by a functional assay, if available.)
7. Grow the culture for 2–3 hours.
8. Remove 5–20 µl of the culture for determination by SDS-PAGE, and harvest the remaining culture by centrifugation. Store the pellets at –70°C.

Note *If induction will be monitored using Coomassie staining, silver staining, or another nonspecific protein staining methods, run a control of CE6-infected BL21 cells harboring the plasmid without a cloned insert.*

TROUBLESHOOTING

Observation	Suggestion
Problems associated with induction time	In certain cases, accumulation of target protein may kill cells at saturation while allowing normal growth in logarithmically growing cultures; in other cases, target protein may continue to accumulate in cells well beyond the recommended 3-hour induction period. To determine the optimal induction period, a time course may be performed in which a small portion of the culture is analyzed by SDS-PAGE at various times following induction
Inclusion bodies	In some cases, protein may form insoluble inclusion bodies at 37°C. In many cases, this protein may be soluble and active if the induction is carried out at 30°C. Inclusion body formation may be used as a purification step by simply spinning out the insoluble material from crude lysates and redissolving the protein in urea or guanidinium-HCl

PREPARATION OF MEDIA AND REAGENTS

Agarose Plates (per Liter) Melt 20 g of agarose in 500 ml of deionized H ₂ O Add the following: 5 g of NaCl 5 g of yeast extract 10 g of tryptone Add deionized H ₂ O to a final volume of 1 liter Autoclave Pour into petri dishes (~25 ml/100-mm plate)	LB Agar (per Liter) 10 g of NaCl 10 g of tryptone 5 g of yeast extract 20 g of agar Adjust the pH to 7.0 with 5 N NaOH Add deionized H ₂ O to a final volume of 1 liter Autoclave Pour into petri dishes (~25 ml/100-mm plate)
LB–Ampicillin Agar (per Liter) 1 liter of LB agar, autoclaved Cool to 55°C Add 10 ml of 10-mg/ml filter-sterilized ampicillin Pour into petri dishes (~25 ml/100-mm plate)	NZY Broth (per Liter) 5 g of NaCl 2 g of MgSO ₄ · 7H ₂ O 5 g of yeast extract 10 g of NZ amine (casein hydrolysate) Adjust the pH to 7.5 with NaOH
NZY Agar (per Liter) 5 g of NaCl 2 g of MgSO ₄ · 7H ₂ O 5 g of yeast extract 10 g of NZ amine (casein hydrolysate) 15 g of agar Add deionized H ₂ O to a final volume of 1 liter Adjust the pH to 7.5 with NaOH Autoclave Pour into petri dishes (~80 ml/150-mm plate)	NZY Top Agar (per Liter) Prepare 1 liter of NZY broth Add 0.7% (w/v) agarose Autoclave NZY Broth with Maltose Prepare 1 liter of NZY broth Autoclave Add 3 ml of a 2 M filter-sterilized maltose solution or 10 ml of 20% (w/v) filter-sterilized maltose prior to use
SOB Medium (per Liter) 20.0 g of tryptone 5.0 g of yeast extract 0.5 g of NaCl Autoclave Add 10 ml of 1 M MgCl ₂ and 10 ml of 1 M MgSO ₄ prior to use Filter sterilize	SOC Medium (per 100 ml) <p>Note <i>This medium should be prepared immediately before use</i></p> 1 ml of a 2 M filter-sterilized glucose solution or 2 ml of 20% (w/v) glucose SOB medium to a final volume of 100 ml Filter sterilize
SM Buffer (per Liter) 5.8 g of NaCl 2.0 g of MgSO ₄ · 7H ₂ O 50.0 ml of 1 M Tris-HCl (pH 7.5) 5.0 ml of 2% (w/v) gelatin Add deionized H ₂ O to a final volume of 1 liter Autoclave	TE Buffer 10 mM Tris-HCl (pH 7.5) 1 mM EDTA

REFERENCES

1. Phillips, T. A., VanBogelen, R. A. and Neidhardt, F. C. (1984) *J Bacteriol* 159(1):283-7.
2. Weiner, M. P., Anderson, C., Jerpseth, B., Wells, S., Johnson-Browne, B. *et al.* (1994) *Strategies* 7(2):41-43.
3. Sambrook, J., Fritsch, E. F. and Maniatis, T. (1989). *Molecular Cloning: A Laboratory Manual*. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY.

MSDS INFORMATION

The Material Safety Data Sheet (MSDS) information for Stratagene products is provided on the web at <http://www.stratagene.com/MSDS/>. Simply enter the catalog number to retrieve any associated MSDS's in a print-ready format. MSDS documents are not included with product shipments.